SERVOX®

unsteril
non-sterile

Tracheo Nose P

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INSTRUCTIONS FOR USE

SERVOX® Tracheo Nose P

Artificial Noses (Heat and Moisture Exchanger HME)

Non-sterile

Symbol explanation

8

Contents

1	Foreword	5
1.1	Intended purpose	5
1.2	Indication	5
1.3	Contraindication	5
2	Safe use	6
3	Description	6
4	Application	6
5	Technical data	7
6	Storage and transport conditions	7
7	Order information	7
8	Disclaimer	8

1 Foreword

EN

This instruction manual is intended for patients, relatives and medical staff. The SERVOX® Tracheo Nose P can be used at home and in a clinical setting. Patients can use the SERVOX® Tracheo Nose P (without assistance from relatives or carers). Application should be in line with the descriptions in this instruction manual and in accordance with the recognised rules of medicine.

Note

Read this instructions for use carefully. They contain important information on the safe and proper use of the product. Keep it for future reference.

A batch number is assigned to the SERVOX® Tracheo Nose P for unique traceability. The batch number is indicated on the packaging next to "LOT".

If you have any questions or problems with the SERVOX® Tracheo Nose P, please contact your field representative or Orbisana Healthcare GmbH, quoting the product's batch number.

1.1 Intended use

The SERVOX® Tracheo Nose P artificial nose is placed on a 15 mm standard connector of e.g. a tracheal cannula or an endotracheal tube. It warms, humidifies and filters the inhaled air by recovering moisture and heat via a filter from the exhaled air, thus reducing irritation and diseases of the respiratory tract.

The SERVOX® Tracheo Nose P must not be used for ventilation.

Note

The SERVOX® Tracheo Noses P is a single-use product and must not be reused or reprocessed, as there is a risk of infection.

1.2 Indication

Use in spontaneously breathing intubated or tracheotomised patients with a tidal volume of 250 – 1000 ml.

1.3 Contraindications

- Severe exsiccosis (dehydration) with excessive secretion of the respiratory tract
- Concomitant use of a speech valve

41294/Version 2022-10-31

2 Safe use

Safe and proper handling is important for your health! Ignoring the safety instructions can have significant effects on your health or even cause life-threatening situations.

Warning: The word "Warning" indicates a situation that, if not avoided, can result in grave injury or even death.

Please observe the following safety instructions:

Warning

Risk of infection: If you use a dirty SERVOX® Tracheo Nose P, you can become infected.

- Replace the SERVOX® Tracheo Nose P immediately if the filter medium becomes heavily contaminated with mucous or dirt, or after 24 hours at the latest!
- Do not use the SERVOX® Tracheo Nose P if it or the packaging is damaged!
- Do not use a fallen SERVOX® Tracheo Nose P, but replace it with a new SERVOX® Tracheo Nose P!
- Do not put filters that have been accidentally ejected (e.g. due to heavy coughing) back into the SERVOX® Tracheo Nose P but use a new SERVOX® Tracheo Nose P!
- Do not use water for rinsing the SERVOX® Tracheo Nose P!
- Dispose of a used SERVOX® Tracheo Nose P immediately after use! **Risk of suffocation/risk of respiratory distress:** If your tracheostoma is blocked accidentally or continuously, you are at risk of respiratory distress or suffocation.
- Never insert the SERVOX® Tracheo Nose P directly into a tracheostoma!
- Do not modify the SERVOX® Tracheo Nose P!
- Immediately remove the SERVOX® Tracheo Nose P if the filter is clogged or in case of increased airway resistance e.g. due to increased accumulation of secretions!
- Never cover or block the filter opening of the SERVOX® Tracheo Nose P with clothing or similar!
- Do not use any atomisers or humidifiers or moisten the SERVOX® Tracheo Nose P with water, as this can increase breathing resistance and block respiration!
- Do not wear the SERVOX® Tracheo Nose P in water, the shower or the bathtub!
- Change the SERVOX® Tracheo Nose P regularly when using in humid environments, such as in the rain!

Increased CO₂ concentration in the respiratory air: A breath volume range above or below 250 – 1000 ml can lead to an increase in CO₂ concentration or insufficient humidification.

- Observe the approved tidal air volume range!


Note

Please inform us and, if necessary, the competent authority if life-threatening illness or injury occurs in connection with this product or if in-patient treatment results from its use.

3 Description

The SERVOX® Tracheo Nose P is a heat and moisture exchanger (HME) that can be connected to a 15 mm standard connector of a tracheal cannula or an endotracheal tube. Moisture and heat of the expired air are preserved in a filter medium upon exhalation and then returned to the breathing air on inhalation. The humidifying effect sets in after just a few breaths. The filter also protects against the ingress of dust, dirt, insects or other foreign bodies with the inhaled air. Using the SERVOX® Tracheo Nose P increases the airway resistance (similar to the physiological resistance of the upper airways), thus supporting the natural lung function.

The housing of the SERVOX® Tracheo Nose P is made from polypropylene (PP) free of latex and DEHP without metallic or electrically conductive materials. There are no restrictions for use during magnetic resonance tomography (MRT). The two filter cores contained in the front of the housing are made of corrugated paper. In the event of a sudden increase in pressure, e.g. due to coughing, the filters are blown out to protect against overpressure in the airways.

SERVOX® Tracheo Nose P

4 Application

Wash the hands thoroughly before use.

Check prior to use: Ensure that the packaging and the SERVOX® Tracheo Nose P are without damage. Do not use if the packaging and the SERVOX® Tracheo Nose P are damaged.

Fitting: Use light pressure to press the SERVOX® Tracheo Nose P with the 15 mm connector onto the standard 15 mm connector of the tracheal cannula or the endotracheal tube until it is firmly fitted.

Removal: Replace the SERVOX® Tracheo Nose after 24 hours at the latest, or immediately in case of heavy mucus or dirt contamination, or if breathing resistance increases. To do this, hold the 15 mm standard connector of the tracheal cannula or endotracheal tube and gently pull the SERVOX® Tracheo Nose P off the 15 mm standard connector and dispose of it immediately in the residual waste.

5 Technical data

The SERVOX® Tracheo Nose P is a single-use medical device.

SERVOX® Tracheo Nose P				
Tidal volume	250 – 1000 ml			
Loss of moisture after 1 hour	10.2 mg/l (for VT= 250ml)			
	12.6 mg/l (for VT= 500ml)			
	17.0 mg/l (for VT= 750ml)			
	17.6 mg/l (for VT= 1000ml)			
Inspiration Pressure drop after 0 hours (VT= 500 ml)	0.43 hPa (for 30 l/min)		Expiration Pressure drop after 0 hours (VT= 500 ml)	
	1.07 hPa (for 60 l/min)			0.32 hPa (for 30 l/min)
	2.13 hPa (for 90 l/min)			0.81 hPa (for 60 l/min)
Inspiration Pressure drop after 24 hours (VT= 500 ml)	0.43 hPa (for 30 l/min)		Expiration Pressure drop after 24 hours (VT= 500 ml)	
	1.11 hPa (for 60 l/min)			0.34 hPa (for 30 l/min)
	2.10 hPa (for 90 l/min)			0.85 hPa (for 60 l/min)
Connections	15 mm standard connector			
Interior volume/dead space	12.5 ml			
Weight	6.8 g			
Maximum period of use	24 hours			

EN

6 Storage and transport conditions

Store or transport in a clean, dry environment, free of dust and protected from light at temperatures between -10°C and 50°C and a relative humidity between 30 %-95 %. Protect from moisture and direct sunlight. Can be used unopened for 3 years from date of manufacture.

7 Order information

Designation	PU	Accessory no.	CPN	Order no.
SERVOX® Tracheo Nose P	30	12.24.07.6055	08717134	28393

41294/Version 2022-10-31

8 Disclaimer

In the event of damage, the manufacturer does not assume liability in case of non-compliance with the instruction manual or the safety instructions or manipulation of the product.

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