SERVOX[®] C15

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INSTRUCTION FOR USE

SERVOX® C15

Closure cap for 15 mm standard connector

EN

Foreword

This instruction for use is intended for patients, relatives and medical staff. Application of the product according to the prescriptions of the attending physician, in line with the descriptions in this instruction manual and in accordance with the recognised rules of medicine! Prior instruction on use by the attending physician or healthcare professional is required. Children or patients with cognitive disorders or with coordination disorders must be supervised and supported during application. For use in home and clinical settings in all age groups.

Each SERVOX® C15 is assigned a batch number (on the packaging next to 'LOT'). Specify these when asking questions about the product.

Note

Read this instruction for use carefully. It contains important information on the safe and proper use of the product. Keep this instruction manual for future reference.

Intended purpose and description

The SERVOX® C15 is used to temporarily occlude a tracheal cannula with a 15 mm standard connector.

When the tracheal cannula opening is closed, the breathing air is directed via the upper trachea and the patient can breathe via Naturalis. Prerequisites for application are a clear upper trachea and sufficient inhalation and exhalation of the patient. Material: Polypropylene (PP). No restriction for use in magnetic resonance imaging (MRI).

Note

The SERVOX® C15 is intended for reuse on a single patient and must not be reused on additional individuals. There is a risk of infection.

Indication

For tracheostomised patients to close the tracheal cannula, e.g. for weaning from the tracheal cannula (decannulation).

Contraindication

- Do not use in laryngectomised patients or on blocked tracheostomy tubes.
- Not to be used in case of allergy to the material polypropylene.
- In case of narrowed upper trachea or granulations in the trachea, use exclusively in combination with screened/fenestrated tracheal cannulas!

Side effects

If the movement of the vocal folds is restricted, a strongly increased inhalation resistance can develop when using the cap, up to an audible stridor (whistling breathing sound). Using the cap may cause an increased respiratory rate and possibly respiratory fatigue or respiratory exhaustion if respiratory performance is impaired.

Safety instructions

Failure to observe the safety instructions may result in serious health effects or even life-threatening injuries.

Warning

Risk of suffocation: Laryngectomised patients and tracheostomised patients with narrowed upper trachea (tracheal stenosis) cannot breathe or breathe poorly when their tracheostoma is closed. There is a risk of respiratory distress or suffocation.

- Use exclusively according to the instructions of the attending physician!
- Never use on laryngectomised patients!
- Never place on blocked tracheal cannulas!
- Use only if the upper trachea is clear! It may be necessary to cough or aspirate beforehand.
- For patients with narrowed upper trachea, use cap exclusively with screened/fenestrated cannulas!
- Before attaching to cannulas that can be blocked or cannulas with exchangeable inner cannula, always check whether the cannula is unblocked or whether a suitable inner cannula is inserted! The cap may obscure the cannula identification when the tracheostomy tube is inserted.
- Always check the patient's breathing and vital parameters during application!
- If breathing suddenly becomes difficult, remove the cap immediately! If breathing difficulties persist, call a doctor immediately!

Warning

Risk of aspiration: Aspiration of objects or food residues may cause respiratory distress or infection.

- Never insert the cap directly into a tracheostoma!
- In patients with increased mucus formation or dysphagia, aspirate any food or fluid residues before using the cap!

Risk of injury: If the cap is deformed, damaged or shows sharp edges, tears or similar, you may injure yourself.

- Do not use even slightly damaged caps, replace them immediately!
- Do not change the cap!

Risk of injury: If you put on or take off the cap with the tracheal cannula already in place, you can injure your trachea by putting too much pressure or pulling on the tracheal cannula!

- Always hold the cannula label and carefully adjust pressure or tension!

Risk of infection: If contaminants or germs get into your tracheostoma, you can get an infection.

- Clean and disinfect the cap before each use!
- Rinse cap thoroughly after each cleaning and disinfection!
- Clean and disinfect caps that have fallen off before putting them back on!

Note

Please inform us or the responsible authorities if you experience a life-threatening illness or injury or permanent physical damage as a result of using our product, or require hospital treatment or have to extend hospital treatment due to using our product.

Application

- Before each use, clean and disinfect according to the instructions for use and check for any damage or deformation.
- To attach, hold the neck flange securely and carefully push the SERVOX® C15 firmly onto the standard connector with light pressure.
- To remove, hold the neck flange and carefully pull off the SERVOX® C15.

Cleaning and disinfection

Attention

- Only ever use suitable and recommended cleaning and disinfection agents! Do not use any aggressive household cleaners or high-proof alcohol! Rinse thoroughly after each cleaning and disinfection!
- Cleaning and disinfection only ever outside of your body!

Required aids and accessories

- 2 plastic cleaning cans with a volume of 400 ml with sieve insert (cannula) cleaning agent for alkali-stable instruments and devices in all medical areas and laboratories (biodegradable surfactants: non-ionic surfactants <5%, anionic surfactants 5-15%, phosphates >30%)
- Cannula disinfectant for plastic tracheal cannulas (e.g. Primediclean (Primed))
- Gloves and access to warm, clean, and running water

Cleaning

- Wash hands.
- Rinse the cap thoroughly with running, lukewarm water.
- Pour a cleaning solution with the cannula cleaning powder according to the manufacturer's instructions in the cleaning tub with the sieve insert. The cleaning solution has to be enough to entirely cover the closure cap during the application time. Wear gloves while handling the cleaning solution in order avoid skin contact with the cleaning solution.
- Place cap in sieve insert and immerse completely in cleaning solution, leaving it to soak for at least 10 minutes (follow manufacturer's instructions).
- Then carefully lift the sieve insert out of the cleaning solution, remove the cap and rinse it thoroughly several times under clear running water (removing any residues of the cleaning solution).
- After cleaning, dry the closure cap with a lint-free cloth or leave to dry completely at room temperature and dispose the cleaning solution in accordance with manufacturer instructions.

Disinfection

- Wash hands.
 - Pour a 2% disinfection solution with the cannula disinfectant according to the manufacturer's instructions in the cleaning tub with the sieve insert. The disinfection solution has to be enough to entirely cover the closure cap during the application time. Wear gloves while handling the disinfection solution in order to avoid skin contact with the disinfection solution. (Observe the manufacturer's instructions).
- Place previously cleaned cap in the sieve insert of the cleaning can, immerse completely and soak for at least 30 minutes (follow manufacturer's instructions).
- Remove the cap and rinse several times and thoroughly under clear running water (removal of any residues of the disinfectant solution).
 - After cleaning, dry the closure cap with a lint-free cloth or leave to dry completely at room temperature and dispose the disinfection solution in accordance with manufacturer instructions.

Storage, transport, service life and disposal

Store or transport disinfected SERVOX® C15 in a clean, dry environment, free from dust and protected from sunlight and light at temperatures between 5°C and 49°C. Multiple use with one patient following the instructions for use for a maximum of 29 days.

Can be stored for a maximum of 3 years in intact packaging if the storage conditions are observed. If no special measures are necessary due to illness (e.g. MRSA), dispose of with household waste.

Disclaimer

In the event of damage, the manufacturer does not assume liability in case of non-compliance with the instruction manual or the safety instructions or manipulation of the product.

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