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INSTRUCTION FOR USE

SERVOCURA®

Hydrophilic Coated Catheter (Ready to use)

Preface

These instructions for use are intended for patients, relatives and healthcare professionals. The product can be used in home care and clinical settings. Consult the attending physician before use! Use the product in accordance with these instructions for use and the approved rules of medicine.

For unique traceability, a batch number is assigned to the product. It can be found on the packaging behind the „LOT“ label. If you have any questions or problems with the product, please contact the sales representative or Orbisana Healthcare GmbH, referring to the batch number.

Note

These instructions for use contain important information for safe and proper use. Keep it for future reference.

Intended use and description

The Hydrophilic Coated Nelaton Catheter – Ready to use is intended to be inserted through the urethra to the bladder for urine drainage, by the patient/user or a health care provider, typically for an individual person who is physiologically incapable of voiding. A Solution for activating coating is included. Do not use water or other lubricant.

Nelaton Tip

Tiemann Tip

Material: Catheter: Thermoplastic polyurethane (TPU), Lubricant: Polyacrylamide derivatives, 2-Propanol, water, Polyurethane

Note

The Product is intended for single use and must not be reused. There is a risk of infection.

Indication

Patients with dysuria, incontinent patients, control of bladder emptying, collection of urine samples or blood clots or sediment from the bladder after surgery.

Contraindication

Acute urethritis, acute prostatitis, acute epididymitis, bleeding in the urinary tract, injury to the urethra, calcareous urolithiasis, allergy against one of the product components.

Side effects / restrictions

Injury or bleeding of the mucous membrane, nosocomial urinary tract infection, shock of the patient, leakage of urine, patients during menstruation have an increased infection risk.

Warnings

- Do not use if the catheter or the individual package has been damaged.
- Do not use beyond expiry date
- Do not re-use. Re-use of the single use product may create a potential harm to the user.
- Dispose of the catheter after use. Reprocessing, washing, disinfection and/or sterilization may compromise product characteristics, causing additional risk of physical harm or infection to the user.

Note

Please inform us and, if necessary, the competent authority if life-threatening illness or injury occurs in connection with this product or if in-patient treatment is required or must be prolonged because of its use.

Use

Note

- Always choose a fitting size
- Do not use if the catheter or its individual package has been damaged
- Catheterize immediately after opening the package
- Remove the catheter when the bladder has been emptied.

Contact your health care professional if you:

- If symptoms such as fever, discomfort when emptying the bladder, frequent urge to urinate occur
- If you notice blood in your urine
- If you feel that the catheter is not emptying the bladder in the way you are used to

Preparation

- Wash your hands thoroughly with soap and water, disinfect them
- Clean and disinfect the area around the urethra following the instructions of attending healthcare professionals
- Check if there is any damage on the product or the packaging
- Spread upper part of package to open. Optional: Use the adhesive tab to attach the package to a dry, clean surface
- Take out the catheter but do not touch the catheter tube
- Pull and adjust the insertion grip to control insertion without having to touch the catheter tube

Insertion of the Catheter

Male

Lift the penis with one hand to straighten out the urethra.

Female

Separate the labia with the thumb and forefinger to locate the urethra. Be sure to distinguish urethra and vagina.

Insert the catheter into the urethra. You can meet resistance at the urethral sphincters. Apply gentle but firm continuous pressure until the catheter advances into the bladder and urine starts to flow.

When urine stops flowing, pull the catheter out about 1-2 cm. If the urine starts flowing again, wait until it stops. Then pull the catheter out again 1-2 cm. When the bladder is empty, remove the catheter and dispose of it according to local regulations.

EN

Storage and transport


Store in a horizontal position in a cool, dry and ventilated place at temperatures between 5°C and 25°C. Avoid contacting with corrosive gases. Protect from rodents and insects.

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Jamjoom Fullcare Corporation
Shanghai International Holding Corp. GmbH (Europe)