HOFFRICHTER OR 5 Oxygen Regulator

HOFFRICHTER
RESPIRATORY CARE

Gebrauchsanleitung
User's Manual

HOFFRICHTER OR 5 Oxygen Regulator

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©2024 Acare Technology Co., Ltd.

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©2024 Acare Technology Co., Ltd.

All rights reserved.

The content published in this user's manual is the sole property of Acare Technology Co., Ltd. All forms of reproduction, editing, distribution and use of any kind, even in part, require the prior written consent of Acare Technology Co., Ltd.

Acare Technology Co., Ltd. reserves the right to amend or replace this user's manual without prior notice.

Please ensure that you are always working with the most current version of this user's manual. Should you have any questions, please contact the provider of the regulator.

Manufactured by:

Acare Technology Co., Ltd.

6F.-3, No.24, Wuquan 2nd Rd., Xinzhuang Dist., New Taipei City 24892,
Taiwan

TEL.: +886-2-2298-8170

FAX: +886-2-2298-8560

Internet: <http://www.acaretech.com>

Declaration of compliance

The system was built in application of the approved quality control system for the design, manufacture and final inspection of the product, and complies with the requirements laid down in Annex II excluding section 4 of the Medical Devices Directive 93/42/EEC and its subsequent amendments (MDD).

With regard to the above, the system bears the CE Marking.

This authorisation is granted by DNV GL PRESAFE SA in its capacity as a Notified Body for CE Marking. DNV GL PRESAFE SA Notified Body identification number is 2460.

Table of contents

Chapter 1	47
Introduction	47
Scope of delivery	47
Symbols	48
Symbols on the device	48
Symbols on the packaging	49
Symbols used in this user's manual	51
Intended purpose / Description of function	52
Contraindications	52
Side effects	52
Chapter 2	53
Safety Information.....	53
Chapter 3	57
Description of Device	57
Chapter 4	59
Commissioning.....	59
Operation	60
Equipment Set-up.....	60
Connecting the regulator to the oxygen cylinder	62
Opening the oxygen cylinder valve.....	63
Adjusting oxygen flow rates.....	65
Closing the oxygen cylinder valve.....	67
Removing the regulator from the oxygen cylinder	68

Chapter 5	69
Cleaning and Reuse.....	69
Chapter 6	71
Maintenance / Inspection.....	71
Chapter 7	73
Service Life and Disposal	73
Chapter 8	75
Annex.....	75
Technical data	76
Warranty.....	77

List of figures

Figure 1: Label packaging.....	49
Figure 2: Description of the regulator	57
Figure 3: Inlet Configuration DIN477 - Oxygen tube at the outlet connection	62
Figure 4: Adjusting knob	66

Chapter 1

Introduction

Scope of delivery

Figure	Name
	HOFFRICHTER OR 5 Oxygen Regulator
	Adapter for connection tube
	User's Manual

Symbols

Symbols on the device

Symbol	Meaning
	CE mark and number from the notified body. The medical device complies with the applicable regulations of EU 93/42/EEC for medical products.
	Medical device: Indicates the item is a medical device.
	Catalogue number: Indicates the manufacturer's catalogue number so that the medical device can be identified.
	Serial number: Indicates the manufacturer's serial number so that a specific medical device can be identified.
	Caution: Indicates that caution is necessary when operating the device or control close to where the symbol is placed, or that the current situation needs operator awareness or operator action in order to avoid undesirable consequences.
	Consult instructions for use: Indicates the need for the user to consult the instructions for use.
	Manufacturer: Indicates the medical device manufacturer.
	Date of manufacture: Indicates the date when the medical device was manufactured.
	Dispose the waste equipment according to the national law.
	Keep away from sunlight: Indicates a medical device that needs protection from light sources.

Symbol	Meaning
--------	---------

Keep dry: Indicates a medical device that needs to be protected from moisture.

Temperature limit: Indicates the temperature limits to which the medical device can be safely exposed.

Symbols on the packaging


Figure 1: Label packaging

Symbol	Meaning
--------	---------

CE mark and number from the notified body. The medical device complies with the applicable regulations of EU 93/42/EEC for medical products.

Medical device: Indicates the item is a medical device.

Catalogue number: Indicates the manufacturer's catalogue number so that the medical device can be identified.

Serial number: Indicates the manufacturer's serial number so that a specific medical device can be identified.

Authorized representative in the European Community/ European Union: Indicates the authorized representative in the European Community/ European Union.

Caution: Indicates that caution is necessary when operating the device or control close to where the symbol is placed, or that the current situation needs operator awareness or operator action in order to avoid undesirable consequences.

Importer: Indicates the entity importing the medical device into the locale.

Consult instructions for use: Indicates the need for the user to consult the instructions for use.

Manufacturer: Indicates the medical device manufacturer.

Date of manufacture: Indicates the date when the medical device was manufactured.

Dispose the waste equipment according to the national law.

Temperature limit: Indicates the temperature limits to which the medical device can be safely exposed.

Keep away from sunlight: Indicates a medical device that needs protection from light sources.

Keep dry: Indicates a medical device that needs to be protected from moisture.

Symbols used in this user's manual

Important information is denoted by symbols in this user's manual. Please ensure that you follow these instructions in order to prevent accidents, personal injury and material damage.

In addition, the local accident prevention regulations and general safety regulations in force in the area of use must be observed.

 DANGER This symbol denotes hazardous situations that lead to serious injuries or death.

 WARNING This symbol denotes hazardous situations that may lead to serious injuries or death.

 CAUTION This symbol denotes hazardous situations that may lead to moderate injuries.

Please note: *Notes provide tips and information for the efficient, correct use of the device.*

Intended purpose / Description of function

The VST Series Click Style Oxygen Regulator reduces pressure for oxygen therapy under the order and direction of a healthcare professional in the clinical environment and in homecare. Instruction of the operator beforehand is necessary for use.

The VST Series Click Style Oxygen Regulator enables pressure-reduced flow controlled extraction from compressed gas cylinders and protects connected devices against overpressure. Setting the flow control knob at the prescribed flow rate adjust the flow rate.

The regulator contains an internal relief valve, which relieves excess pressure, if a regulator malfunction occurs. The relief valve operates as a safety mechanism; it is not intended to maintain the downstream pressure or flow within normal limits. In the event of the malfunction, excess pressure will be vented either intermittently or continuously by this safety feature. Never close or cover the outlet DISS 9/16-18UNF, centered on the regulator.

If the situation occurs, close the cylinder valves immediately and remove the regulator from service until repaired.

Contraindications

There are no known contraindications for the use of the HOFFRICHTER OR 5 Oxygen Regulator.

Side effects

Possible side effects from the use of the HOFFRICHTER OR 5 Oxygen Regulator are not known.

Chapter 2

Safety Information

Please read this instruction of use thoroughly before operating equipment and follow directions carefully.

 DANGER **Danger of fire and explosion!**

- Make sure that the person installing or using the regulator has clean hands that are free of hand lotions, petroleum based product, hair spray or similar products. In combination with fat, oxygen can ignite itself.
-
- Inspect the cylinder thoroughly for oil and grease. If oil or grease is visible or suspected, do not use the cylinder. Do not attempt to clean or repair the cylinder. Inform the oxygen supplier of this condition immediately.
 - Inspect the regulator thoroughly for oil and grease. If oil or grease is visible or suspected, do not use the regulator. Do not attempt to clean or repair the regulator.
 - Do not expose oxygen or oxygen equipment, as the regulator to fire, heat, sparks, electrical switches or other possible sources of ignition.
 - Do not smoke near oxygen or oxygen equipment.
 - Do not store or use oxygen equipment at temperatures of less than -18°C or more than 60°C and near heat sources generally.
 - Do not cover the cylinder under clothing or in any way, which could trap oxygen and so cause a fire.
 - If you detect any leakage of the oxygen cylinder, immediately turn handle clockwise until the cylinder valve is closed. If leakage persists, place the cylinder outdoors and inform the gas supplier immediately.
 - Not closing the cylinder valve before removing the regulator will cause a rapid depressurization of the cylinder. This can cause a fire.
-

 **WARNING** Danger of injury!

- Always stand behind the cylinder outlet or regulator when opening the cylinder valve.
- Always open the cylinder valve slowly and pressurize the regulator gradually.
- Turn the cylinder valve fully off before removing the regulator.
- The device cannot be operated by blind and severely visually impaired persons.

Please note: *We ask you to inform us and, if needed, the competent authority if you suffer a life-threatening illness or injury or permanent bodily harm in connection with the use of our product or if you need to be treated as an in-patient or prolong an in-patient treatment due to the use of our product.*

Chapter 3

Description of Device

The following figure explains the individual components of the pressure reducer.


Figure 2: Description of the regulator

A Adjusting knob B Display flow rate C HOFFRICHTER OR 5 Oxygen Regulator
D Screw connection for adapter to the connecting tube E Pressure display
F Connection for oxygen cylinder with connecting nut

Chapter 4

Commissioning

Operation

Please note: *It is very important to keep product in original packaging for 12-24 hours to acclimatize to room temperature before use.*

The operating and storage temperature for the regulator should reflect typical environmental conditions of a medical facility environment. Do not store or use oxygen equipment at temperatures of less than -18°C or more than 60°C.

Equipment Set-up

Follow these steps to set up a filled cylinder:

1. Always support and fix the cylinder to secure it from falling or rolling over in a place with adequate ventilation.
2. Clean hands and make sure, that your hands are free from hand lotions, petroleum based product, hair spray or similar products.
3. Make sure, the valve outlet of the cylinder is always facing away from you or any person during operation.
4. Remove the protective dust cover (if present) from the cylinder valve. Keep the cover for protection of the cylinder valve during storage or transportation, in case, the regulator must be removed.
5. Inspect the cylinder valve for traces of dirt, dust, oil or grease, especially on the face of the cylinder valve where the regulator connects. Do not use the cylinder, if oil or grease is visible or suspected. Dirt or dust must be removed safely with a lint-free-cloth. Be sure there is no debris left on the cleaned surface.
6. To purge any undetectable debris from the valve seat, open the oxygen cylinder valve slightly, then re-close the valve.

Follow these steps to set up the regulator:

-
1. Remove the protective dust cover from the inlet nut.

 2. Inspect the regulator for traces of dirt, dust, oil or grease, especially on the connection. Do not use the regulator, if oil or grease is visible or suspected.

 3. Screw the included connector for the oxygen supply tubing onto the regulator outlet connection.

Connecting the regulator to the oxygen cylinder

1. Screw the regulator inlet connection nut onto the cylinder valve clockwise.
 2. Tighten the nut tightly, do not use any tools to do so.
 3. Attach the oxygen supply tubing to the regulator outlet connection.
-


Figure 3: Inlet Configuration
DIN477 - Oxygen tube at the
outlet connection

A Adapter for tube
B Oxygen tube

Opening the oxygen cylinder valve

WARNING Risk of injury!

- Always stand behind the cylinder outlet or regulator when opening the cylinder valve.
 - Always open the cylinder valve slowly and pressurize the regulator gradually.
 - If the regulator, cylinder or other devices, connected do not operate normally, stop operating, immediately.
 - Never try to remove the regulator before the cylinder valve is closed and all pressure has been relieved.
-

1. Stand behind the cylinder valve so that the cylinder is between you and the regulator.
 2. Before opening the oxygen cylinder valve, make sure, that the oxygen supply tubing between the regulator and the patient is not longer than 210 cm and not obstructed kinked or blocked, this can cause backpressure.
 3. Turn the flow control regulator in the off position, before opening the oxygen cylinder valve. This can be accomplished by turning the adjusting knob counterclockwise to the "0" flow position (in the flow display). Do not force the knob beyond this point.
 4. SLOWLY turn the cylinder valve handle counterclockwise to open the valve. The valve should be open approximately one turn (i.e. 360°).
-

-
5. If a hissing sound occurs, there is the leak in the cylinder valve, regulator to cylinder connection, or in the regulator. In this case, close the cylinder valve and follow the following directions.
- ⇒ Close the cylinder valve and open the regulator shortly to relieve all pressure by turning the adjusting knob clockwise and wait for the oxygen flow to stop. Then return the adjusting knob on the regulator again counterclockwise to the "0" flow position.
 - ⇒ Slowly turn the cylinder valve handle counterclockwise to open the valve. The valve should be open approximately one turn (i.e. 360°).
-

 CAUTION Risk of injury due to defective cylinder valve!

Should the hissing sound continue, close the cylinder valve and remove the regulator.

- ⇒ The sealing washer or the regulator is leaking and must be removed from service and returned for repair.
 - ⇒ The cylinder valve is leaking. Place the cylinder outdoors and inform the gas supplier immediately. Do not attempt to use the cylinder.
-

Adjusting oxygen flow rates

CAUTION Danger due to excessive back pressure!


Backpressure in excess of 2 psi will cause inaccuracy of the oxygen flows that are indicated on the adjusting knob. Always make sure, that the attached supply tubing is not too long and not kinked

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Do not cover or block the vent hole in regulator body with labels, decals, or other adhesive material. Blocking the vent may inadvertently affect regulator preset pressure, resulting in inaccurate flow rates.

Adjusting the adjusting knob on the end of the regulator determines the oxygen flow of the regulator.

Please note: At back pressure the flow rates are inaccurate. The adjusting knob must be set at a marked flow setting. Placing the knob between settings will result in inaccurate flow. Although the adjusting knob "clicks" into the specified flows on the knob, oxygen does flow out of the regulator in between the "clicks" of the knob.

A Adjusting knob

Figure 4: Adjusting knob

For oxygen therapy applications, turn the adjusting knob clockwise (flow) until the flow on the knob is indicating the flow prescribed by your attending physician.

Closing the oxygen cylinder valve

1. Always close the oxygen cylinder after use.
 2. With the oxygen still flowing, turn the cylinder valve or wrench clockwise until it is closed.
 3. Allow the gas pressure in the regulator to escape, completely. (Wait until the manometer shows no pressure.)
 4. After all pressure has been relieved, turn the regulator to the off position by turning the adjusting knob counterclockwise to the "0" flow position. Do not force the knob beyond this point.
-

Removing the regulator from the oxygen cylinder

WARNING Danger due to unlocked cylinder valve!

Always close the cylinder valve completely before removing the regulator.

Please note: A cylinder showing less than 300 psi on the high pressure gauge should be exchanged for a full cylinder to ensure an adequate supply of oxygen.

It is not necessary to remove the regulator unless the cylinder is being exchanged for a new one.

If a humidifier cup is used, disconnect it from the regulator before removing the regulator from the cylinder to prevent water entering the regulator.

-
1. Be sure the cylinder valve is closed.
 2. Open the regulator shortly to relieve all pressure before removing it from the cylinder.
 3. Lose the nut by turning it in a counterclockwise direction. Do not use tools.
 4. Once loose, turn the nut until it is disengaged from the cylinder valve. Be sure to support the regulator during removal.
-

Chapter 5

Cleaning and Reuse

⚠ CAUTION In order to use the device for another person, please follow the instructions specified in the section Cleaning.

⚠ CAUTION **Danger due to incorrect cleaning agents!**
Do not use aggressive or corrosive cleaning solutions.

⚠ CAUTION **Danger due to water in the regulator!**
Never let water enter the regulator or submerge the regulator in any form of liquid.

1. Disconnect all connections before cleaning.
 2. Clean the exterior of the oxygen regulator after each removal from the oxygen cylinder, with a clean and lint free cloth, dampened with mild detergent and water.
 3. Afterwards, wipe it dry with a clean lint free cloth.
 4. Store the oxygen regulator in a clean area free from grease, oil and other sources of contamination or reuse as like it is described in "Chapter 4" on page 59.
-

Chapter 6

Maintenance / Inspection

The oxygen regulator should be checked and repaired when malfunction occurs or is suspected. Every 5 years the oxygen regulator must undergo a periodic security inspection. Security inspections and repairs must be performed by a qualified personnel. Please contact the customer service of Hoffrichter GmbH. Never attempt to repair any part of the regulator or cylinder valve by yourself.

Chapter 7

Service Life and Disposal

The expected lifetime of the regulator is 10 years. Dispose the oxygen regulator in accordance with the local regulators. Please contact the customer service of Hoffrichter GmbH or a certified scrap recycler.

 **WARNING Biohazard!**

Potential for biohazard. Clean the regulator before disposing.

Chapter 8

Annex

Technical data

The VST 485 Series Click Style Oxygen Regulator is suitable for regulating the pressure of standing oxygen cylinders with a maximal pressure 200 bar.

Weight	434,7 g
Dimensions	500 mm x 223 mm
Suitable for	Oxygen (O ₂)
Oxygen supply tubing	210 cm, Ø = 4,9 mm
Operating pressure	4,5 bar
Outlet pressure nominal	4,5 bar +/-0,34 bar
Selectable outlet flow	1–5 l/min, adjustable in the steps 0 – 0,5 – 0,8 – 1 – 1,5 – 2 – 2,5 – 3 – 4 – 5 l/min
Accuracy	+/-20 %
Inlet	High pressure manual connection → G3/4" thread DIN 13260
Inlet configuration	DIN 477
Inlet pressure range	200 bar
Outlet connection	9/16-18 UNF
Serial number	See engraving (8-digit device number) The first 4 numbers date year and month of manufacturing
Manometer	In bar, scaling is in steps of 10 bar

Operation Environmental limits	-18 °C – 60 °C
Storage Environmental limits	-20 °C – 60 °C
Humidity	Max. 95 %, non-condensing
Maintenance	Every 5 years and in case of malfunction / suspected malfunction.

Warranty

This unit is warranted under normal use against any and all manufacturing defects from date of purchase for a period of 1 year.

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HOFFRICHTER GmbH
Surpass Biomedical Europe B.V.